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NECK MASSAGER

Developed based on the principles of traditional Chinese medicine. Modern FAR infrared technology penetrating even into deeper tissue layers combined with warmth, magnetism, electric impulses and vibration massage. Ergonomic design for best wearing comfort.

PRODUCT SPECIFICATIONS

Electric impulse: 2 automatic + 4 manual therapy modes

Vibration massage: 6 modes

Magnets

Operated by remote control with an illuminated LCD display

PLEASE NOTE!

Important! Before first use load up completely (approx. 4 hours)!

Operation: Connect the adapter to the massager, then connect the massager to the power supply. The device will be fully charged in about 3-4 hours.

Application:

1. Apply the device to clean skin only. Slightly moisten the back of your neck and the metal with a cloth.

2. Put the power switch into the "ON" position and place the device onto the back of your neck.

3. Press the  button to start the operation. You will hear a beep when the device and the remote control have been synchronised (**as shown in the picture, the infrared element and the device need to be in contact with each other [SYNCHRONISATION!]!**). By pressing the  button again, the application will be stopped.

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 Increase

5. Selection of the different modes

By operating the buttons provided in the  selection mode. 6 different modes can be selected.

1. Automatic mode: encompasses vibration + warmth; warmth + electric impulse; warmth; vibration; electric impulse; automatic circulation

2. Vibration and warmth

3. Warmth and electric impulse

4. Warmth

5. Vibration

6. Electric impulse

6. Adjusting the intensity

Selection of the  symbol: Control via the buttons (section 4)

7. Adjusting the warmth

Selection of the  symbol: Control via the buttons (section 4)

8. Adjusting the electric impulses

Selection of the  symbol: Control via the buttons (section 4)

9. Adjusting the vibration

Selection of the  symbol: Control via the buttons (section 4)

10. Timer function

Selection of the  symbol: Control via the buttons (section 4)

11. Battery

Connect the device to the power supply via the adapter, and switch the ON/CHARGE button to CHARGE. The POWER display will flash in green. As soon as it stops flashing and turns to permanent light (after approx. 3-4 hours), the device will be fully charged. Then all functions can be used for approx 5 days (for the recommended two applications daily for 15 minutes each). If the device is used more often, the battery life is shortened.

Suggestion for the application:

Twice daily 15 minutes each, with vibrations and electric impulses as well as 15 minutes warmth.

CAUTION

1. Do not use this device on persons with a pacemaker.
2. Do not use it close to microwaves or other infrared-emitting devices.
3. Pregnant women and persons with cardiovascular diseases should consult a doctor before using the device.
4. Do not use during your sleep or under the shower.
5. May not be used by children or persons who cannot communicate.
6. Not suited for persons with sensitive skin or heat-sensitive persons.
7. An increase of pressure, output, or warmth will not lead to a positive therapeutic effect.
8. Do not disassemble the device, plunge it into water, or clean it by directly applying water.

CLEANING AND CARE:

1. Use a moist cloth for cleaning.
2. Store in a dry, well ventilated place.
3. Keep out of direct sunlight and do not store in an extremely moist or hot place.

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PLEASE NOTE:

In case the remote control does not work though the main unit has been charged and both units of the device are in their correct position, please exchange the coin batteries (CR2032) of the remote control.

TECHNICAL SPECIFICATIONS

Name: Neck massager

Model no.: KTR-101

Weight: 630 g

Power: Device: 6VDC - 6 W; remote control: CR2032 coin battery 3VDC

Electric impulse frequency: 1-1000 Hz

Electric impulse range: 20-400 µs

Heating temperature: +40° ~ +48°C (with an ambient temperature of +25°C)

DISPOSAL: Do not dispose this product as unsorted municipal waste, collection of such waste separately for special treatment is necessary.

Important: Please make sure the battery is disposed of properly, complying with current legislation. Note on disposal: Discharged batteries do not belong in household waste. Some batteries may be recycled, and may be accepted for disposal at your local recycling centre.

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