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Spannung: 230V ~50Hz

Leistung: 22W (100W max.)

BACK AND NECK HEATING PAD
Instruction Manual

IMPORTANT SAFETY INSTRUCTIONS RETAIN FOR FUTURE USE

Read the instructions

Do not use when folded or creased

Do not insert pins

Not to be used by very young children (0-3 years)

Keep carefully in case it is needed!

This electric heating pad is a high-quality product designed to last for many years if it is properly used and cared for. Read the following instructions carefully to prevent damage or danger from incorrect use. Keep the instructions so that you can refer back to individual points when necessary.

If you let somebody else use the pad, you must be sure to let it have the instruction on use.

This appliance is not intended for medical use in hospitals.

- Please read carefully and follow the instructions.
- Use the device only for intended use.
- Only connect to the voltage shown on the rating plate.
- The appliance must not be used when folded or crumpled.
- This appliance must not be used by persons insensitive to heat and other very vulnerable persons who are unable to react to overheating.
- Do not use the heating pad for persons who cannot move without help, small children, persons who cannot feel heat, or animals. Persons who can feel heat to some extent may only use the pad with supervision.
- Do not operate the heating pad unsupervised.
- Caution! Never fall asleep while the heating pad is in operation.
- Prolonged use can cause burns on the skin.
- Do not pierce the appliance with needles, pins or other sharp objects.
- Do not tug or twist the power lead.
- Do not expose the switch or power lead to moisture or liquids.
- The manufacturer or similar qualified persons must replace if the cord is damaged, in order to avoid a hazard.
- Do not use the appliance when it is damp or wet.
- The device may not be switched on using a time switch.
- The switch should not be placed on or beneath the appliance or be covered when the appliance is in use.
- Avoid sharp bending or bundling of the appliance.
- While the appliance is switched on:
 - Do not place any large or heavy objects on it such as suitcases, etc.,
 - Do not place any source of heat on it such as hot water bottles, heating pads or suchlike.
 - Do not allow the switch to be concealed under blankets, pillows or the like and lie on the appliance.
 - Do not fold the pad when hot.
- Examine the appliance frequently for signs of wear or damage. If there are such sign or if the appliance has been misused or does not work, return it to the supplier before switching on again.
- Children should be supervised to ensure that they do not play with the appliance.
- This appliance is not to be used by young children over the age of three unless the controls have been pre-set by a parent or guardian or the child has been adequately instructed on how to operate the controls safely.
- Children below the age of 3 years are not allowed to use this device due to their inability to react to overheating.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- When the appliance is likely to be used for an extended period (falling asleep), the control should be adjusted to setting 1 for continuous use.
- The switch or control unit must not be allowed to get wet during washing and during drying the cord must be positioned to ensure that water does not flow into the switch or control.
- When the appliance is not in use, the plug should be removed from the mains and the appliance should be stored dry, without additional weight on it.
- When not in use, store as follows:
 - When storing the appliance, allow it to cool down before folding.
 - Do not crease the appliance by placing items on top of it during storage.
- Please use only with the plug of the type 102C3.

SWITCHING ON

Please always ensure that the heating pad is not bundled or folded when it is in use. The heating pad has been specially designed for use on back and neck. You can adjust the elastic fastener in the lower part of the heating pad to fit the pad to your body.

On initial use the heating pad may smell of synthetic material. However, this odor will disappear after a short while.

The appliance may only be used with 230 Volts a.c. Switch on the pad using the hand-operated switch. When switch on it have 3 times beep –beep –beep sounds to ensure the user that the switch has been operated correctly.

4 = maximum warmth
3-2 = moderate warmth
1 = minimum warmth
0 = appliance off

3 COLORS LED

The heating pad provides 3 colors LED for using friendly.

Orange – Heating up

Green – Keep warm

Red and Beep sound – The red color and beep sound give warning for the below situation, when the RED led lights on please switch off the pad immediately.

- When the power is on and detachable connector is unplugged. The beep sound of interval for 1 sec. Please connect the detachable cable into the pad then switch on the controller again.
- When critical components are damaged: the beep sound will be continuously. Please disconnect the plug from wall outlet immediately.
- When the inside of wiring is broken or other components failure, it cannot work anymore. The beep sound of interval for 1 sec.
- When the heating pad is folded and full covered, the pad will reach to higher temperature. The beep sound of interval for 2 sec. Please switch it off and check if incorrect using.

Above mentioned warning cannot cover all situation of defect as some of components defect cannot give red light or beep sound.

QUICK HEATING

The heating pad with a quick heating function that user can feel warm within 7-10 minutes when turn on the switch at high setting.

TIMER FUNCTION

Now you can experience the pleasant warmth. After about 180 minutes the pad silently switches off automatically. In order to put the heating pad into operation again you must first set the switch to the '0' position. After about 5 seconds you can switch on the heating pad again.

Due to safety concern, the back and neck heating pad will reduce temperature to setting 1 automatically after using 90 min.

CAUTION:

The electronic components inside the switch may lead to a slight warming of the switch when the pad is on. Therefore the switch must not be covered or placed on top of the heating pad when it is in operation.

You should not rely on using your hands to gauge the pad temperature because they are much less sensitive to heat

than the surface of the human body. A temperature, which may feel tepid to your hand, may feel unpleasantly hot on your body.

Qualified persons using original replacement parts may only undertake repairs to the device. Repairs carried out by unqualified persons may cause substantial risk to the user.

CLEANING AND CARE

Please follow below symbols as indicated on the washing label.

- Do not dry clean
- Do not bleach
- 30° C extra gentle washing program
- Do not tumble
- Do not iron

HEATING PAD

- Before cleaning the pad, unplug it from the mains outlet and let it cool down.
- The heating pad is equipped with a detachable cable. Disconnect the plug and remove the cable from the heating pad.
- Please note that the heating pad should not be dry-cleaned, wrung out, machine-dried, pressed or ironed.
- Set washing machine to a gentle washing cycle at 30°C. The heating pad has to be washed only together with other textiles. Use a washing powder or liquid for delicate fabrics and dose it according to the manufacturer's instructions.
- Please note that excessive washing will wear out the heating pad. For this reason, the heating pad should only be washed a maximum of 10 times throughout its entire life.
- Only use the appliance when it is complete dry.

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