Pro Anti-Cellulite Vacuum Device InnovaGoods® - Wellness Beauty

FRANÇAIS

immerse the appliance or cable in liquid. Do not handle it with wet hands Do not use near water or damp areas. Suitable for domestic use only. Use the device for its intended purpose only. Store the device in a dry, dust-free place. Do not leave the appliance unsupervised when it is turned on. If the appliance, the cable or the pieces looked damaged or do not work correctly, do not use them or attempt to repair them; they should only be repaired by a qualified technician. This device is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by an adult responsible for their safety. This device is not household waste. Please follow the environmental regulations in place for its recycling

ENGLISH

DEUTSCH

ESPAÑOL

ITALIANO

CONTENTS

- 1 ergonomic Pro Anti-Cellulite Vacuum Device InnovaGoods® 800 mAh
- 1 head with double roller
- 1 head with single roller
- 1 power charger DC 5 V
- 1 USB Cable

INSTRUCTIONS

- Plug in the Pro Anti-Cellulite Vacuum Device InnovaGoods®.
- Slide the on/off switch to the ON position and press the START button.
- Pro Anti-Cellulite Vacuum Device InnovaGoods® comes with two heads:
 - A. The head with one roller is designed to be used on arms, stomach and thighs.
 - B. The head with two rollers is designed to be used on legs, waist and glutes.
- Insert the head of your choice into the device placing the ends together so that they meet, and then turn it clockwise.
- Use the side button to adjust the intensity of the Pro Anti-Cellulite Vacuum Device InnovaGoods®.
- Before each treatment, apply reducing oil or body oil on the skin to enhance the massage and avoid any discomfort.
- For best results, massage each area for 5 to 8 minutes.
- To change movement or roller head, press the switch at the back to stop the suction.
- It can be used without a cable for areas that are difficult to reach, approximate duration: 8 minutes.

WARNINGS

If you feel any discomfort when using the device, discontinue use and consult a doctor. After using the appliance you may notice redness on the treated skin; this is normal and it will disappear quickly. Do not use the device on areas where the skin is infected or sore. To avoid the build up of hard-to-remove residue, clean the device after each use. Fully charge the battery before its first use. Make sure that the voltage of the power source corresponds to that of the appliance. Do not pull on the cable and ensure it does not tangle. Keep the cable away from hot and sharp surfaces. Clean with a damp cloth. Do not use abrasive products for cleaning. Do not

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