English

Congratulations on your purchase of Velform® Vacu Clear, the pore cleaning device that will give you glowing skin. Velform® Vacu Clear gently eliminates facial blackheads and clogged pores without squeezing. Its powerful suction mechanism removes dirt and oil from the skin.

CONTENTS

- 1- Warnings
- 2- Product Components
- 3- Battery Installation
- 4- How to Change Suction Heads
- 5- Instructions for Use
- 6- Cleaning and Maintenance

1. WARNINGS

SAFETY INSTRUCTIONS

READ AND FOLLOW ALL WARNINGS AND INSTRUCTIONS BEFORE USING THIS PRODUCT. INJURY CAN RESULT FROM IMPROPER USE. KEEP THESE INSTRUCTIONS IN A SAFE PLACE.

- Do not wear necklaces, earrings or any facial jewellery while using this product. Necklaces, earrings and other types of jewellery should be removed before using Velform® Vacu Clear.
- This product is not for use on areas of the body that have undergone any sort of surgery or have any scars.
- If you experience pain while using this product, discontinue use immediately.
- This product is not waterproof and should be not used in the shower or any situation that involves water.
- Do not store this product in the bathroom or a room that has high moisture content.
- Do not operate in the presence of explosives and/or other flammable fumes or liquids.
- Do not store or use this product near sources of heat or fire.
- Keep out of reach of children.
- Keep batteries away from children (harmful if swallowed).

2. PRODUCT COMPONENTS

ON/OFF and Low/High Selection LOW/HIGH (Figure 1).
Suction head (Figure 2).
Battery compartment (Figure 3).

Small head (Figure 4):
Gentle suction; perfect for everyday use. Good for use on sensitive skin and nose area.

Large head (Figure 5):
Strong suction; perfect for blackheads and deep pore cleaning.
Note: Do not use around the eyes.

Oval head (Figure 6):
Oval opening for getting into corners around the nose and mouth, and for sensitive skin around the eyes.

Sonic head (Figure 7):
Helps exfoliate dead skin while pulling dirt and oil from the skin.

3. BATTERY INSTALLATION

1. Open the battery compartment cover (Figure A).
2. Insert 2 AA batteries (not included) into the battery compartment. Make sure to line up the (+) and (-) markings on the battery with the (+) and (-) markings on the inside of the battery compartment (Figure B).
3. Reinstall the battery compartment cover.

4. HOW TO CHANGE SUCTION HEADS

1. Place one hand on the suction head and the other hand on the unit (Figure C).
2. Firmly pull off the suction head and reinstall your desired suction head in the same manner (Figure D).

5. INSTRUCTIONS FOR USE

1. Start by washing your face with warm water, steam or a hot wet towel to open up your pores.
2. Make sure to dry your face and hands completely before moving on to the next step.
3. Select your desired suction head (see point 4 on how to change suction heads). Velform® Vacu Clear can be turned on by sliding the ON/OFF switch up. The middle position is low suction and the top position is high suction (Figure F).
4. Apply the Velform® Vacu Clear so the head comes in contact with your face (Figure G). Slowly move the Velform® Vacu Clear around your face without leaving it in one spot for more than two seconds. For best results, focus on your problem areas.
5. Note: DO NOT move the device back and forth on a specific area.

Helpful Tips

- After using Velform® Vacu Clear, make sure to apply moisturizer to the treated area.
- Before and after use, do not use facial scrubs or other exfoliating products.
- Do not press down during use, given that the Velform® Vacu Clear's suction capabilities are enough to clean your pores efficiently.
- Slight redness is normal and will disappear after approximately 30 minutes.

Velform® Vacu Clear can be used once a week for dry skin. For oily and mixed skin, it can be used twice a week. Do not use it for longer than five minutes a day.

6. CLEANING AND MAINTENANCE

After each use, clean the suction head by rinsing it with warm water. Allow the suction head to air dry. Do not attempt to clean the main unit by rinsing it with water because it is not water-resistant. To clean the main unit, wipe it down with a soft, damp cloth.

CAUTION: DISPOSAL OF MATERIALS

A symbol of a crossed-out wheeled bin means you should find out about and follow local regulations about disposing of this kind of product.
Do not dispose of this product as you would other household waste.
Dispose of this device in accordance with the corresponding local regulations.
Electrical and electronic devices contain hazardous substances that can have harmful effects on the environment and/or human health and should be recycled properly.

Australia & New Zealand only

This product has the benefit of certain Consumer Guarantees. These are prescribed by the Australian Consumer Law & the New Zealand Consumer Guarantees Act 1993 both of which provide protection for consumers. There is no express warranty for this product in Australia or New Zealand. The above paragraph refers to other countries.

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English

Congratulations on the purchase of Velform® Vacu Clear, the pore cleanser for a radiant skin. Velform® Vacu Clear gently removes blemishes and dirt from the pores on the face, without the need for pressure. Its powerful suction mechanism removes dirt and oil from the skin.

- 1- Warnings and safety instructions
- 2- Product parts
- 3- Inserting the batteries
- 4- Replacing the suction heads
- 5- Instructions for use
- 6- Cleaning and maintenance

- 1- WARNINGS AND SAFETY INSTRUCTIONS. READ AND FOLLOW ALL WARNINGS AND INSTRUCTIONS BEFORE YOU USE THIS PRODUCT. IMPROPER USE CAN CAUSE INJURY. KEEP THESE INSTRUCTIONS in a safe place.
- Do not wear necklaces, earrings or other jewelry while using the device. Remove necklaces, earrings and all other jewelry before using Velform® Vacu Clear.
- This product should not be used on areas of the body that have had surgery or have scars.
- If you experience pain while using the product, stop immediately.
- This product is not waterproof and should not be used in the shower or other areas with water.
- Do not use the product in the bathroom or other humid areas.
- Do not use the product near heat sources or open flames.
- Keep the device out of children's reach.
- Keep batteries out of children's reach (choking hazard)

- 2- PRODUCT PARTS
- Intensity selector: HIGH/LOW (Figure 1).
- Suction tip (Figure 2).
- Battery compartment (Figure 3).

- Small tip (Figure 4):
- Light suction, ideal for daily use. Suitable for sensitive skin and the nose zone.

- Large tip (Figure 5):
- Strong suction, perfect for comedones and for deep cleaning of pores.
- Note: Do not use around the eyes.

- Oval tip (Figure 6):
- Oval opening to reach the nostrils and the mouth corners, as well as for the sensitive area around the eyes.

- Sonic tip (Figure 7):
- Helps exfoliate the skin and remove dirt and oil.

- 3- BATTERY INSERTION
- 1. Open the battery compartment cover (Figure A).
- 2. Insert two AA batteries (not included) into the battery compartment. Make sure the polarity symbols (+) and (-) on the battery match the symbols (+) and (-) inside the compartment (Figure B).
- 3. Replace the cover.

- 4- REPLACE THE SUCTION HEADS

- 1. Hold one hand on the suction tip and the other on the device (Figure C).
- 2. Pull firmly on the tip to remove it and insert the other tip in the same way (Figure D).

- 5- INSTRUCTIONS FOR USE
- 1. Wash your face with warm water, so that the pores are better open; keep your face over steam or moisten it with a warm, damp cloth.
- 2. Dry your face and hands well before continuing with the next step.
- 3. Choose the suction tip you want to use (see paragraph 4 - Replacing the suction tips). Switch the Velform® Vacu Clear on by sliding the power switch up.
- 4. Place the Velform® Vacu Clear so that the tip touches your face (Figure G). Move the Velform® Vacu Clear slowly over your face and do not keep it in one place for more than two seconds.
- 5. Note: Do not move the device forward or backward in the same direction.

- Practical tips
- After using Velform® Vacu Clear, apply a moisturizing cream to the treated area.
- Do not use facial peels or any exfoliating products before or after using the device.
- It is not necessary to press the device against the skin, as the suction of Velform® Vacu Clear is sufficient to clean the pores effectively.
- Small red spots may appear, but they will disappear within a few minutes.

- We recommend that you use Velform® Vacu Clear once a week on dry skin. Use the device twice a week on oily or combination skin. Do not use the device for more than 5 minutes per day.

- 6- CLEANING AND MAINTENANCE
- After each use, clean the suction tip by rinsing it with warm water. Let the tip dry in the air. Do not clean the main unit under running water, as it is not waterproof. Simply use a damp cloth to clean the device carefully.